

# True Clarity

## APEX® Multifocal Corneal GP Lens

The Apex Multifocal offers a modern approach to presbyopic contact lens fitting. Designed with a unique back surface platform to improve comfort, centration and response. The lens also features a completely different way of delivering the near, intermediate and distance optics. And, it's backed by our comprehensive fit warranty.

### The APEX Multifocal Parameters


**Base Curve:** Custom

**Diameter:** Any (default: 10.5mm)

**Power:** Custom

**Add:** Custom Constructed for the Pupil Size

• Available in a Back or Bi toric



Optimal Fit

### A Modern Approach

- Empirically ordered or trial sets are available.
- Large diameter for improved centration and stability
- High aspheric back surface to improve alignment and comfort.
- Unique multifocal optics to provide quality near, intermediate and distance vision
- Created specifically for the individual's power and pupil size.
- Easy to follow fitting system supported by our consultants

### Diagnostic Selection

BASE CURVE	DIAMETER	POWER	ENGRAVING
7.10	10.5	-6.00	A
7.20	10.5	-5.50	B
7.30	10.5	-5.00	C
7.40	10.5	-4.50	D
7.50	10.5	-4.00	E
7.60	10.5	-3.50	F
7.70	10.5	-3.00	G
7.80	10.5	-2.50	H
7.90	10.5	-2.00	I
8.00	10.5	-1.50	J
8.10	10.5	-1.00	K

**+2.00 Add @ 3 mm Pupil – Diagnostic Rings @ 3 mm and 5 mm**

\* Contact your Account Manager for more information or pricing.

Your Partners in Visual Excellence.

**pts** PRECISION  
TECHNOLOGY SERVICES

**Cardinal**  
Contact Lens Inc.

#### Western Canada

**Precision Technology Services Ltd**

866 Cordova Street East, Vancouver, BC,

Canada V6A 1M4

Trevor Kojima: **604 215 4206**

Toll Free: **1 800 663 4248**

Email: [trevor@ptsoptics.com](mailto:trevor@ptsoptics.com)

#### Eastern Canada

**Cardinal Contact Lens Inc.**

#7-38 Cowansview Rd., Cambridge, Ontario,

Canada N1R 7N3

Amy Chow: **437 993 2020**

Toll Free: **1 800 263 3973**

E-mail: [amy@cardinalcontactlens.com](mailto:amy@cardinalcontactlens.com)